

# Esthetic Procedure for Photon Facial

By Lynn Dellapasta  
Westshore Skin and Health Center  
Westlake, Ohio



- 1) Cleanse face thoroughly to remove oil and makeup using gentle cleanser. Rinse well and pat dry.
- 2) Turn off bright overhead lights to allow WARP 10® as much brightness as possible with no competing light.
- 3) Using 2 WARP 10s (one in each hand) start on the neck and hold each in place on either side of the throat. Deliver two 80 second doses to the neck in that area. Wait for the delay light to go out after the first dose, then simply press the start button again for the second dose-this is very easy.
- 4) Move to the next area on the neck moving upward towards the jawline always delivering 2 doses in each area being treated.
- 5) Keep treating the lower face (along the jawline too) with 2 doses per area.
- 6) Pay attention to where you are treating because it is not necessary to overlap areas.

7) Continue moving your WARP 10s in an upward pattern and delivering 2 doses per area. When you reach the forehead and finish the entire surface of the forehead, you have completed your full treatment.

8) Finish the treatment by either applying a quality mask (with antioxidant ingredients) and leave in place for 7 - 8 minutes or simply apply an anti-aging serum w/an oil-free sunscreen afterwards. If you choose to use a mask, rinse off well and then apply a serum and sunscreen afterwards.

9) Procedure is complete at this point.

We recommend doing this 2-3 times per week for 3 weeks in a row. After that schedule is finished, have client come in once or twice a month for maintenance.